

COVID-19 QUICK DEPLOYMENT BUNDLE

Environmental Monitoring

AeroScout® RTLS and Hugs® (MobileView®) users can add Environmental Monitoring, including tags and remote installation support.



Challenges

As hospitals respond to the increased patient load resulting from the COVID-19 pandemic, there has been a surge in vital materials that must be stored within prescribed temperature ranges, including test kits, pharmaceuticals, and blood products. Storage conditions must be monitored continually to ensure patient safety and protect these valuable materials.

Solution

STANLEY Healthcare's AeroScout® Environmental Monitoring solution is a wireless hospital temperature monitoring solution that helps ensure the safe storage of pharmaceuticals, vaccines, laboratory samples and blood products.

Easy to deploy with no fixed wiring, the AeroScout solution leverages your Wi-Fi network to cover every part of the hospital. Our state-of-the-art temperature sensors provide local and remote alerting for temperature excursion events, temperature monitoring to -200°C and integrated contact closure to enable notifications of door ajar events. Data are sent to a centralized RTLS platform for later analysis and regulatory compliance. Our sensors log data in case network connectivity is lost, and release data for upload when the connection is re-established.

How It Works



TAG

Sensor tags are placed in refrigerators or freezers



MONITOR

Temperature is monitored over Wi-Fi



ALERT & CORRECT

Real-time alerts when condition levels deviate, and easily log corrective actions




REPORTING

Historical data for regulatory compliance and operational reporting and analytics

COVID-19 RESPONSE LIMITED-TIME OFFER*

To support the needs of our customers as they work to reduce the spread of COVID-19 while caring for infected patients, we have developed a special time-limited offer that enables you to quickly deploy our wireless Environmental Monitoring solution.

This Environmental Monitoring bundle is available to current AeroScout RTLS and Hugs (MobileView) customers.



Deployment

We have designed a simplified and remote implementation process that enables hospitals and healthcare facilities to quickly take advantage of our RTLS Environmental Monitoring solution. This enables you to be fully operational on an accelerated schedule.

Conditions

- Offering is intended for existing STANLEY Healthcare customers without Environmental Monitoring but with **MobileView 5.4** already installed
- Installation location is limited to Pharmacy, Blood Banks, and Labs
- Deployment is to cover refrigerators and freezers only; no ambient monitoring
- Temperature tags are configured and installed by the customer. An Apple iOS device is required to configure tags
- Support services, including customer familiarization with the solution, will be delivered remotely
- Requires remote access by STANLEY Healthcare to the MobileView server
- Includes Software Maintenance for the first year
- Out of Range alerts are limited to a single layer, with no escalation
- Device categories are limited
- Glycol for the temperature probes to be provided by the customer
- Tags are configured in unidirectional mode

CONTACT US

Contact your STANLEY Healthcare Client Executive to learn more or call
+1-888-622-6992

INTERNATIONAL

+44 (0)1753 500608

*DISCLAIMER

All product and pricing offerings are subject to a mutually-agreed quote and acceptance by STANLEY Healthcare of customer's purchase order. STANLEY Healthcare reserves the right to set criteria for this promotion based on factors including, but not limited to, contractual terms and product availability. The offer period begins on April 1, 2020 and expires on June 30, 2020. Customer must submit a purchase order to STANLEY Healthcare no later than June 30, 2020. This offer may not be applied retroactively for orders placed prior to April 1, 2020. STANLEY Healthcare reserves the right to modify or cancel this offer at any time without notice. Offer valid for customers only in the United States, Canada and other selected geographies.



About STANLEY Healthcare

Over 15,000 hospitals and senior living communities rely on STANLEY Healthcare solutions to empower caregivers to deliver better care. STANLEY Healthcare is a part of Stanley Black & Decker and a proud supporter of the Alzheimer's Association®. Learn more at stanleyhealthcare.com.